



UC Davis MMPC-Live Protocol

Feces Collection for Gut Microbiome Analysis

Version: 1.0

Revision Date: 10/20/2023

Replaces version:

Edited by: Michael Goodson & Trina Knotts - UC Davis Metabolism & Metabolic Health Core

[Summary](#)[Reagents and Materials](#)[Protocol](#)[Reagent Preparation](#)

Summary:

The work of Gordon and colleagues (Nature. 2006 Dec 21;444(7122):1027-31) has shown that obesity can result in marked shifts in the gut microbiome in mice and other models including humans. While the role of the microbiome remains to be fully elucidated, the gut microbiota can no longer be ignored as a potentially important factor when assessing metabolic phenotype. This service involves sequencing the **16S gene variable region (V3-V4 using primer pair 515F-806R) pair end sequencing (2x300bp) by Illumina (20,000 avg seq read depth) of feces, cecal, or other GI contents**. Microbial genomic DNA extraction, library preparation and Illumina MiSeq PE300 sequencing of samples collected by the Metabolism and Metabolic Health core at the UC Davis MMPC-Live center are performed by the UC Davis Host Microbe Systems Biology Core (<https://health.ucdavis.edu/medmicro/hmsbcore/index.html>). Sequences can be processed through a bioinformatics pipeline (Qiime2) to taxonomically classify them and to assess alpha and beta diversity of the community. In addition, Principal Components Analysis (PCA) or partial least squares- discriminant analysis (PLS-DA) can leverage variances in the relative microbial abundances to better understand how specific microbes contribute to separation by group. Correlational analyses can identify which variables of host metadata associate with specific microbes. The Core's gut microbiome assay will employ this approach to uncover unique microbiota fingerprints in test mice. One caveat is that with fecal samples, patterns are only a surrogate for actual gut microbiota patterns, and may not exactly reflect the intestinal populations.

Reagents and Materials:

<i>Reagent/Material</i>	<i>Vendor</i>	<i>Stock Number</i>
Tweezers		
Disinfectant solution (70% isopropanol or 70% ethanol)		
Labeled pre-weighed 1.5 ml tubes		
Gloves		
Small weigh boats		
Measuring cup/beaker for food weight measurements		
Large beaker for mouse weight measurements		

Protocol:

1. Mice are singly housed in clean cages with fresh food and water and low bedding.

IMPORTANT: Special consideration should be given to the type of bedding used during the 48 hour collection period. When bedding is used, we recommend Carefresh Bedding (Absorption Corp, #868744), using slightly less than the normal amount of bedding per cage (low bedding conditions). For the short collection period, no bedding or special raised cage bottoms will simplify the fecal collection process; however, these modifications will require IACUC amendment approval. Also, corncob bedding (i.e. Bed- O'Cobs, 1/8in; Andersons Lab 8B) is not advised due to the difficulty in sifting through it to collect feces and the fact that it has caloric value (which can complicate analyses- the mice may eat it).

2. Weigh mouse and "food-in" (grams of rodent diet provided at start of experiment).
3. Save several pellets of rodent diet for analysis (e.g., for energy or macronutrient content)-Store at 4°C in airtight container.
4. At 24h, weigh mouse and "food-out" (grams of rodent diet remaining)
5. At the end of 48h collection period, weigh mouse and "food-out".
6. Swap out the cage bottoms with new clean bottoms, replace feeder insert/water bottle and return animals to housing system.
7. For each cage, carefully inspect and remove each piece of bedding using tweezers to collect all feces. Inspect each pellet of feces to assure that there is no spilled food or fibers of hair or bedding attached. If significant food spillage is found in cage bottom, weigh in pre-tared weigh boat to add back to food measurement.
8. Place feces in labeled pre-weighed 1.5 ml tube.
9. Weigh tubes and record weight.
10. Subtract this value from original tube weight to determine fecal "wet" weight.
11. After collection, fecal samples are stored at -80°C until analysis.

Note: It is not critical to perform food intake measurements during this collection time but it is beneficial for analysis/QC.

Reagent Preparation:

None.