



UC Davis MMPC-Live Protocol

Preparation of Fecal Extracts for ELISAs

Version: 1.0

Revision Date: 10/20/2023

Replaces version: None

Edited by: Michael Goodson - UC Davis Metabolism & Metabolic Health Core

[Summary](#)

[Reagents and Materials](#)

[Protocol](#)

[Reagent Preparation](#)

Summary:

This protocol is used to create a dispersion of feces or other gut contents in extract buffer to allow for the release of soluble analytes like lipocalin or calprotectin for subsequent analysis by ELISA or other assays.

Reagents and Materials:

Reagent/Material	Vendor	Stock Number
Walter Stern (100E) 4mm Glass Beads	Fisher Scientific	11-312B
BioSpec Zirconia/Silica Beads (one pound bottle) 0.1mm diameter	Fisher Scientific	NC0362415
2 ml screwcap tubes		
1x PB		
MPBio Fastprep-24 5G homogenizer, or equivalent	MP Biomedicals	116005500
Microcentrifuge capable of 12,000 x g		

Protocol:

1. Prepare a 2 ml screwcap tube with ~3 mm deep 0.1mm dia. silica beads and one 4 mm dia. glass bead. This is roughly equivalent to the MPBio lysing matrix E, which would clearly also work for this assay.
2. Weigh prepared homogenization tube from step one.
3. Add 1-2 mouse fecal pellets (~8-50 mg wet weight) to a prepared homogenization tube.
4. Weigh homogenization tube with fecal material to determine fecal weight.
5. Add 500 μ l Extract Buffer to each homogenization tube/feces sample.
6. Secure tubes and homogenize for 40 s at 6.0 M/s at room temperature in the MPBio Fastprep-24 5G homogenizer (human feces protocol).
7. Remove tubes and centrifuge for 5 min at $>12,000\times g$ at 4°C.
8. Transfer supernatant to a fresh microfuge tube and assay. Store remaining homogenate at -80°C.

Reagent Preparation:

- **Extract Buffer:** 1x PBS + 0.1% Tween-20