



# Alkaline Phosphatase Activity

Version: 1

Edited by: Gary Cline, John Stack: Yale MMPC Analytical Core

**Summary:** Procedure to measure the amount of Alkaline Phosphatase activity. Alkaline Phosphatase (ALP) activity is measured from the hydrolysis of 4-nitrophenylphosphate to 4-nitrophenoxide ion (monitored at 405 nm) and phosphate.

## Reagents and Materials:

Reagent/Material	Vendor	Stock Number
Alkaline Phosphatase Reagent	Prolabs	R85120
Assayed Control Serum 1	Prolabs	R83082
Assayed Control Serum 2	Prolabs	R83083

**Protocol:** Analysis by automated system Cobas Mira Plus.

- 1) Calibrate Cobas for Alkaline Phosphatase Activity analysis by running two assayed control serum.
- 2) Sample handling as performed by the Cobas Mira Plus.
  - a) Pipette 3  $\mu$ L of sample into a cuvette slot.
  - b) Add 150  $\mu$ L of Alkaline Phosphatase Reagent.
  - c) Mixture is incubated at 37°C and spun for 10 minutes.
  - d) Absorbance is measured at 405 nm.

## Reagent Preparation:

**Alkaline Phosphatase Reagent:** Add the appropriate amount of water (6.5mL) to the reagent bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

**Assayed Control Serum 1:** Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

**Assayed Control Serum 2:** Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.