



Thiobarbituric acid reactive substances (TBARS) Assay

Version: 1
Replaced by version: N/A
Edited by: Eva Feldman, M.D.

[Summary](#)
[Reagents and Materials](#)
[Protocol](#)
[Reagent Preparation](#)

Summary: Plasma concentrations of thiobarbituric acid reactive substances (TBARS) are an index of lipid peroxidation and oxidative stress. The protocol describes how the AMDCC quantitates TBARS in the animal models.

Reagents and Materials:

Reagent/Material	Quantity Required	Vendor	Stock Number
Thiobarbituric Acid (TBA)	ICN		190284
Trichloroacetic Acid	Sigma		490-10
Malonaldehyde bis(dimethyl acetal)	Fisher		AC 14861-1000

Protocol:

Reagent Preparation:

Thiobarbituric Acid (TBA): 67mg in 1mL DMSO then add 9mL H₂O.

10% Trichloroacetic Acid (w/v): in H₂O.

1,1,3,3-tetramethoxypropane: 4.167 μ L in 1mL Ethanol then add 49mL H₂O. (500 μ M)

Sample Preparation:

Plasma:

- Place 100 μ L plasma into a labeled 1.5mL micro-centrifuge tube.

Tissue:

- Label 1 sets of 1.5mL micro-centrifuge tubes, 1 set screw top tubes and 1 set of 0.5mL tubes.
 - Weighed out ~20mg and sonicate in 200 μ L RIPA buffer + inhibitors.
 - Sonicate.
 - Centrifuged @ 3000 for 10 min @ 4⁰.
- Remove 10 μ L aliquot into the 0.5mL tubes for protein analysis.
 - Place 100 μ L lysate into a labeled 1.5mL micro-centrifuge tube.
 - Add 200 μ L ice cold 10% Trichloroacetic acid to precipitate protein.
 - Incubate for 15 minutes on ice.
 - Prepare standards as follows:

CONCENTRATION (μ M)	H ₂ O	TETRAMETHOXYPROPANE
0	500	-----
0.625	500	500 from tube 3
1.25	500	500 from tube 4
2.5	500	500 from tube 5
5.	500	500 from tube 6
10	800	200 from tube 7
50	500	500 from tube 8
100	800	200 of 500uM stock

- Centrifuge samples @ 2200 x g for 15 min. at @ 4⁰C.
- Place 200 μ L supernatant and standards into new labeled screw top 1.5ml tube.
- Add and equal volume of 0.67% (w/v) TBA.
- Incubate in a boiling water bath for 10 min.
- Cool. Sample is ready for assay.

Performing Assay:

- While samples are cooling, layout on computer and save as TBxxxxxx.sed where xxxxxx is the date in yyddmm format.
- Load 150 μ L into each standard well in duplicate.
- Load 150 μ L into each samples well in duplicate.

4. Put in plate reader and press start.

Reading the Plate:

- Record absorbance at 532nm.
 1. Turn on Multiskan and open your saved file TBxxxxxx.sed.
 2. Place plate onto Multiskan holder and click
 3. Select Process>Curve Fit. Choose the appropriate data (usually Measure1), then click .
 4. Save Curve Fit data sheet as an Excel file into the Data folder/TBARS data folder. Use the naming convention TBxxxxxx.xls, where xxxxxx is the date in yymmdd format.