



## *UC Davis MMPC-Live Protocol*

### **Non-Esterified “free” Fatty Acids (NEFA)**

Version: 1.0  
 Revision Date: 9/17/2012  
 Replaces version: none  
 Edited by: Peter Havel - UC Davis Metabolism & Metabolic Health Core

[Summary](#)  
[Reagents and Materials](#)  
[Protocol](#)  
[Reagent Preparation](#)

### **Summary:**

The Wako enzymatic method relies upon the acylation of coenzyme A (CoA) by the fatty acids in the presence of added acyl-CoA synthetase (ACS). The acyl-CoA thus produced is oxidized by added acyl-CoA oxidase (ACOD) with generation of hydrogen peroxide, in the presence of peroxidase (POD) permits the oxidative condensation of 3-methy-N-ethyl-N( $\beta$ -hydroxyethyl)-aniline (MEFA) with 4-aminoantipyrine to form a purple colored adduct which can be measured colorimetrically at 550 nm.

### **Reagents and Materials:**

<i>Reagent/Material</i>	<i>Vendor</i>	<i>Stock Number</i>
Calibrator	Wako	276-76491
Reagents	Wako	999-34691 995-34791 991-34891 993-35191
Microplate		
Platereader		

### **Protocol:**

1. Reconstitute Color Reagent A with 50 ml of Solvent A and Color Reagent B with Solvent B.
2. Add 5  $\mu$ l of calibrator and sample to each well.
3. Add 200  $\mu$ l of Reagent A to each well. Incubate at 37°C for 5 minutes. Read at 560 nm.

**IMPORTANT:** Make sure not to add any bubbles to the wells when dispensing reagents, this will interfere with reading in the platereader.

4. Add 100  $\mu$ l of Reagent B to each well. Incubate at 37°C for 5 minutes. Read at 560 nm.
5. Subtract blank readings from final readings. The assay will be linear so the unknown samples can be calculated as (Sample Absorbance  $\div$  Calibrator Absorbance)  $\times$  Calibrator Concentration.

## **Reagent Preparation:**

**Reagent A** – reconstitute Color Reagent A with Solvent A

**Reagent B** – reconstitute Color Reagent B with Solvent B